



Simplicity in every analysis.
Power in every insight.

POWERFUL | EASY TO USE | RELIABLE | COMPLIANCE-READY



Waters™

Your go-to platform for
biopharmaceutical product
attribute and process
analysis.



1

Introduction

2

waters_connect
Ecosystem

3

Usability and
Reliability

4

Proteins and
Peptides

5

Nucleic Acids
and
Delivery Vehicles

6

Bioprocess
Monitoring

7

Configure Your
System

The BioAccord LC-MS System: Where **Simplicity** meets **Power**

Simplicity of operation combined with a purposefully designed hardware and informatics platform provides users of all levels of LC-MS experience the **Power** to understand their product and process and confidently measure their control.



Reduce system downtime and improve separation quality with compatible MaxPeak™ Premier Systems and Columns.



Align with your sustainability goals with environmental impact metrics. Look for labels in 2025.



Achieve data integrity and regulatory compliance goals with a compliance-ready platform and our Global Services Team.



Enable new and experienced users with support from the Waters Global Services team to deploy the BioAccord System and support biopharmaceutical workflows.



Integrated Informatics Workflows from Characterization to QC Release

The waters_connect™ informatics ecosystem provides the flexibility to deploy LC-MS systems for product and process analysis, attribute monitoring, and product release. Compliance-ready biopharmaceutical applications streamline and automate data acquisition, processing, review, and reporting.

Characterization

Attribute Monitoring

Routine/QC Release

Attribute Characterization LC-MS platforms

Xevo™ MRT

Xevo™ G3 QTof MS



Users/Client PC(s)



Network Server
Managed Data Scientific Library
Compliance Tools



Attribute Analysis LC-MS systems

BioAccord LC-MS System



LC Optical System(s)
and Client PC



Peptide
MAM



UNIFI™



INTACT
Mass



MAP
Sequence




CONFIRM
Sequence




Bioprocess
Monitor

Designed for Usability and Reliability

LC-MS complexity has traditionally restricted the technology to central labs and experienced users. The BioAccord LC-MS System has been designed to provide robust operation and communicate system readiness to produce high quality biopharmaceutical data.



Commercial calibration and system check standards support automated system setup and regular system verification checks to ensure proper operation. Detailed feedback on the front panel guides a user to quickly resolve detected issues.



BioAccord System robustness was evaluated using ~2000 intact IP-RP-LCMS analyses of a ~7.8 kDa PS 25-mer oligonucleotide (GEM91) over the course of 120 days of system operation. Data analysis was automated using the waters_connect INTACT Mass Application, and instrument setup-calibration (Cross hair Circles) was conducted on Days 1, 35 and 75. Average mass error of 7.6 ppm from the calculated mass was observed, and 3 sigma control limits (99.7% of data) were calculated at 0.12 Da.




"The ease of use and system robustness enables less-expert users to perform LC-MS tasks, providing greater productivity for advanced analytical strategies such as Multi-Attribute Monitoring (MAM)"

Arnaud Delobel, R&D and Innovation Director, Quality Assistance


[WATCH TESTIMONIAL](#)

The BioAccord System for Protein-based Therapeutics and Vaccines



Accessible attribute analysis for protein & peptide therapeutics, from discovery screening to commercial manufacture.




- Protein Identity and Glycovariants
- Bioconjugation Ratios
- Denatured and Native analysis




- LC-MS^E Peptide Mapping for Attribute Definition
- Peptide MAM for Attribute Monitoring



- Subunit Identity and Variation
- Localization of Modifications
- Glycovariant Profiles



- Released N-glycan profiling using RapiFluor-MSTM, 2AA, 2-AB, or your preferred label



"The BioAccord System has accelerated our native analysis of proteins, increasing throughput threefold. An increase in mass range expands our ability to investigate AOCs, mAbs aggregates, multi-specific mAbs, and oligomer therapeutic proteins, enabling us to support increasingly complex therapeutics"

Dr. Oscar Hernandez Alba,
BioOrganic Mass Spectrometry
Laboratory (LSMBO), IPHC, CNRS,
Strasbourg University



Complete workflows from sample prep to automated data analysis



Unparalleled application-specific consumables, automation, and system-check standards for optimal performance and reproducibility



Purposely designed compliance-ready informatics with intuitive application workflows

waters_connect[™]



UNIFI[™]



INTACT
Mass




Peptide
MAM

The BioAccord System for Nucleic acid-based therapeutics and vaccines


Streamlined attribute analysis and impurity analysis for therapeutic and vaccine development and commercialization.




- Intact Mass Confirmation
- MS2 Sequence Confirmation
- Impurity Identification and Quantification



- LNP Component and Impurity Profiling
- Impurity Detection and Characterization




- Sequence Confirmation
- mRNA 5' Capping Efficiency
- mRNA 3' polyA Tail Heterogeneity



- Confirmation of Viral Particle Product Identity
- Profile and Map Subunit Modifications

siRNA cRNA
ASO saRNA
gRNA LNP
mRNA AAV
tRNAⁿ



"An ease of operation which has clearly never been seen before in other kinds of instrumentation"

*Davy Guillaume, Research Associate,
University of Geneva*

[WATCH TESTIMONIAL](#)




Complementary workflows for holistic analysis of nucleic acids and delivery vehicles



Unparalleled application-specific consumables, automation, and system-check standards for optimal performance and reproducibility



MAXPEAK.
PREMIER



Purposely designed informatics with intuitive application workflows

waters_connect™



UNIFI™



INTACT
Mass



MAP
Sequence



CONFIRM
Sequence


The BioAccord System for Bioprocess Analytics

Simplify in-process analytics, monitor CQAs and CPPs and gain critical process insights to accelerate process development.




Out of the box connectivity for sample-to-result automation

- Ambr® 250 High Throughput and Ambr® 15 system connectivity
- Bioprocess Walk-Up Solutions



Collate data, improve DOE and generate critical process insights



Monitor in-process CQAs and CPPs on a single platform

- Identify CQAs
- Measure abundance
- Monitor trends with time course data

Titer
Glycans
(Intact, Reduced, Subunit)
Purity
(Aggregates/ fragments)
Peptide Monitoring
Cell Culture/
Spent Media



"The synergy between the BioAccord LC-MS System and the Ambr® bioreactor significantly enhances the efficiency of upstream process"

Dr. Li Wang, Director of BioDev Analytical Sciences, WuXi Biologics



Improve operational efficiency and analytical robustness with a readily-deployable platform solution



SARTORIUS

Unparalleled application-specific consumables, automation, and meaningful informatics connections



waters_connect™

Purposely designed informatics with intuitive application workflows



INTACT Mass



Bioprocess Monitor

Flexible: Configure your BioAccord System for your analytical needs

The BioAccord System is deployable from discovery to QC/Product Release and used for a wide variety of biopharmaceutical product and process analytical workflows. Configure your systems for the analyses you employ and the throughput you require.

Optical Detection

TUV/PDA

Add an FLR for Labeled Glycan Analysis

Sample Management

Variable or Fixed Loop Sample Injection



Configurable option to add Automation Portal Sample Interface



Configurable option to add Sample Organizer for Higher Throughput Analysis (up to 7680 samples)

Solvent Management

Binary UPLC (BSM) for Ultimate Gradient Fidelity

Quaternary UPLC (QSM) for Greater Mobile Phase Flexibility



Column Management

Column Manager or Column Heater



Configurable option to add two Auxiliary Column Manager(s) for up to four additional columns



Configurable option to add 30 cm Column Heater/Cooler for LC columns up to 30 cm length and up to 7.8 mm diameter

Vacuum Pump System

Configurable option for Oil-based or Oil-Free (Air-Cooled) Vacuum Pump System for reduced Maintenance and Improved Environmental Sustainability



Made for Biopharma Analysis

- Developed and Tested with Biomolecules
- Installed and Trained with Biomolecules



Your trusted partner for success

- Outcome-Based Training with Dedicated Application Kits
- Highest customer satisfaction score for service and support among all instrument vendors.*



waters.com/bioaccord

Visit videos.waters.com to see the
BioAccord LC-MS System in action

For your local sales
office, please visit
waters.com/contact



Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
waters.com

Waters™

Waters, BioAccord, waters_connect, MaxPeak, Xevo, UNIFI, UPLC, and *RapiFluor*-MS™ are trademarks of Waters Technologies Corporation. All other trademarks are the property of their respective owners.

Ambr is a registered trademark of Sartorius Stedim Biotech GmbH or its affiliates. JMP is a registered trademark of SAS Institute, Inc

©2025 Waters Corporation. May 25-14005 720007339EN